




Animas® Vibe™ Features Guide

In comparison to the OneTouch® Ping®* insulin pump.

Feature	Animas® Vibe™
Battery Insertion/ Verify Screen	<ul style="list-style-type: none"> Time, mode and date must be edited in order to confirm and move past the Verify Screen.
Main Menu	<ul style="list-style-type: none"> The main menu includes CGM† option.
Contrast Button/ CGM† Shortcut	<ul style="list-style-type: none"> When using the Dexcom G4® PLATINUM Sensor and Transmitter with the pump, pressing the  button while pump is in sleep mode will awaken the pump to one of the CGM Trend graphs or the CGM Data screen. While the pump is "awake", pressing the  button will adjust the contrast of the display screen.
Max Prime	<ul style="list-style-type: none"> There is a 20 unit maximum prime amount. You will be unable to prime more than 20 units at one time. After rewinding, pump displays Disconnect Set before the other instructions as an additional safety notification.
ezCarb and ezBG Bolus	<ul style="list-style-type: none"> On the Bolus Total Screen, one press of the  populates the blinking yellow cursor field with the suggested bolus amount.
Food List	<ul style="list-style-type: none"> The Food List is listed on the ezCarb Home screen, and provides an easy and accurate way to obtain carb totals when using the Bolus Calculator.
CGM† Data Screen/CGM† Trend Graphs	<p>View CGM data on your pump using the CGM Data and CGM Trend screens:</p> <ul style="list-style-type: none"> CGM Data screen - provides a snapshot of the current glucose status and Insulin on Board (IOB). CGM Trend screens - provides a more detailed view of the glucose status, and includes a graphical display of the CGM readings over a selected time period (1, 3, 6, 12, or 24 hours).
CGM† Alerts	<p>Set the pump to display and sound (beep/vibrate) an alert when:</p> <ul style="list-style-type: none"> CGM readings are outside the target range (High Alert / Low Alert). CGM readings may be rising or falling too quickly (Rise Rate / Fall Rate). The Transmitter is not within RF range of the pump. There are other CGM problems that require attention.
CGM† History	<p>Historical CGM records can be reviewed on your pump, including:</p> <ul style="list-style-type: none"> Start date and time of the current CGM session. Date, time and value of the last BG calibration entered in the pump. Date, time, alert code, and description of at least 300 of the most recent CGM alerts.
Advanced SETUP Screen Changes	<ul style="list-style-type: none"> The I:C, ISF, and BG Target are on the SETUP Adv 1 screen. The Audio Bolus setup is on the SETUP Adv 7 screen.
Status Menu	<ul style="list-style-type: none"> The Status menu includes a Pump Information screen (Status 6) that displays the pump serial number and CGM transmitter ID/serial number. Status 7 displays additional pump codes (previously on Status 6)
Navigation Arrows	<ul style="list-style-type: none"> Navigation arrows are located at the bottom of the History, SETUP, and Status menus to allow forward and backward movement through the menus.
Other	<ul style="list-style-type: none"> If the pump buttons are locked (Tamper Resistant Feature), the pump will display a screen with additional instructions for unlocking the buttons. Meter-Remote included on OneTouch® Ping® model only.

*This comparison of features applies only to the OneTouch® Ping® insulin pump. For information about features of other Animas pumps please visit www.animas.com.

†Continuous glucose monitoring (CGM) with Animas® Vibe™ is indicated for adults 18 years and older only.

420-465-00 Rev A, ANM-14-4136C, © 2014 Animas Corporation 2014/08 Rx Only.

